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INTRODUCTION


DIGI-FLEX® is a unique hand and finger exercise system designed for you to work each finger individually on its own spring loaded button.

DIGI-FLEX® develops isolated finger strength, flexibility and coordination as it builds hand and forearm strength. There are five (5) models of graduated resistance. You will gain a positive sense of accomplishment as you advance from color to color.

The exercises shown in this pamphlet are illustrative only. Your doctor or therapist should approve all rehabilitative exercises and recommend the number of repetitions and intensity of use. Notify them if any pain is experienced.

This pamphlet has been created in association with:

Joanne Kassimir, OTR, CHT
LONG ISLAND HAND REHABILITATION
275 Hillside Avenue
Weston Park, NY 11598

DIGI-FLEX® SYSTEM resistance scale:

	Per Finger (A)	Entire Unit (B)
Yellow	15 lbs	50 lbs
Red	30 lbs	100 lbs
Green	50 lbs	160 lbs
Blue	70 lbs	230 lbs
Black	90 lbs	310 lbs

(A) Pressure required to fully depress a single button to the center bar.

(B) Pressure required to fully compress the entire unit (all seven (7) springs).

ATTENTION THERAPIST—All pound values are accurate to plus/minus 4 oz. (.25 lbs).

HOOK POSITION

Place fingertips on individual buttons and position ergonomic base with the hook and over the web space between the thumb and index finger. Keeping fingertips slightly bent, flex all fingers toward the center of the unit.

GROSS GRASP

Place the middle of the fingers on the buttons with the fingertips wrapped over the top. Positioning the hook of the ergonomic base over the web space between the thumb and index finger, flex the thumb and all fingers toward the center of the unit as if making a fist.

TIP TO TIP PINCH

Place the tip of the thumb on the center of the ergonomic base. Position the finger to be exercised on one of the center buttons. Keeping both tips slightly bent, flex the thumb and the finger being exercised toward the center of the unit at the same time.

KEY PINCH

Place the hook of the ergonomic base on the inside of the bent index finger. With the tip of the thumb slightly bent, flex the thumb on the first button.

Note: For increased resistance, flex the thumb tip on two (2) buttons at the same time.

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
TRIPOD PINCH

Place the thumb tip in the center of the ergonomic base with the index and long fingertips resting on the two center buttons. Keeping the fingertips bent, flex the thumb and fingers toward the center of the unit at the same time.

—B—

FINGER EXTENSION

Place the unit on the table with the buttons facing you. Rest your forearm on the tabletop and place the face of the nail of the finger you wish to exercise against the tip of the button. Holding the unit in place with the unaffected hand, extend (straighten) the finger pushing the button with the face of the fingernail.

—B—

TRIGGER PINCH

Loosely hold the unit vertically with the buttons toward the fingers and the hook of the ergonomic base positioned over the web space between the thumb and index finger. With the middle of the index finger, press the second button as if squeezing a trigger.

Note: This can also be done using the above position and the tip of the index finger.

—B—

FINGERTIP FLEXION

Place all fingertips on their own button with the tip of the thumb on the center of the ergonomic base. Keeping all fingertips bent, flex the thumb and all the fingers toward the center of the unit at the same time.

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INTRINSIC PLUS POSITION

Place ergonomic base at the back of the thumb with the buttons resting at the base of the fingers. Hold the tip of the thumb and fingers straight out. Flex the thumb and all fingers toward the center of the unit at the same time.

—12—

THUMB FLEXION

THUMB TIP FLEXION: Place the middle of the fingers on the ergonomic base and the tip of the thumb on the second or third button (whichever is more comfortable). Flex the thumb down toward the fingers.

THUMB BASE FLEXION: Follow the above instructions except place the base of the thumb on the second or third button.

Note: For increased resistance, place thumb in the center of the second and third buttons and flex.

—13—

THUMB OPPOSITION


Place end of ergonomic base at the bottom of the small finger with the base of the thumb resting on the end button. Flex thumb in a downward motion toward the bottom of the small finger.

Note: To keep the unit steady, hold the other end with the unaffected hand.

—14—

THUMB ADDUCTION

Place unit flat on the table. Position hand so the base of the thumb rests on the first button and one end of the bar rests against the side of the base of the index finger. While holding the unit in place with the unaffected hand, slide the thumb toward the index finger to squeeze the unit.

Note: The hand must remain flat on the table for maximum effectiveness. Be sure not to use the thumb tip to do this motion.

—15—